

Supporting Information

for

Multicomponent reaction access to complex quinolines via oxidation of the Povarov adducts

Esther Vicente-García¹, Rosario Ramón¹, Sara Preciado¹ and Rodolfo Lavilla^{*1,2}

Address: ¹Barcelona Science Park, Baldiri Reixac 10–12, 08028, Barcelona, Spain and ²Laboratory of Organic Chemistry, Faculty of Pharmacy, University of Barcelona, Avda. Joan XXIII sn, 08028, Barcelona, Spain

Email: Rodolfo Lavilla - rlavilla@pcb.ub.es


*Corresponding author

Experimental details


Table of contents

Synthetic procedures and characterization data for new compounds.....	S2
¹ H and ¹³ C NMR spectra	S4
Particle size analyses of MnO ₂ samples	S6


9-Bromo-5-(4-chlorophenyl)-3,4-dihydro-2H-pyrano[3,2-c]quinoline (18)



3-(6-Bromo-2-(4-chlorophenyl)quinolin-3-yl)propan-1-ol (19)



3-[2-(4-Chlorophenyl)-6-methylquinolin-3-yl]propanamide (22a)




MnO₂ Particle size study:

All the manganese dioxide samples were analyzed with a LSTM 13 320 series Laser diffraction particle size analyzer to determine the particle size.


Reagent code	d ₁₀ (μm)	d ₅₀ (μm)	d ₉₀ (μm)	<75%	<95%	<10 μm (%)
310700	1.475	4.302	11.16	7.022	14.53	87.5
243442	77.66	138.4	173.2	156.8	183.2	1.57
217646	0.577	4.240	16.04	9.583	19.69	75.3
213490010	1.509	7.555	37.30	18.16	51.39	57.3
138-09675	2.743	25.70	63.10	44.77	73.05	21.6
Old Sample	0.932	11.46	47.78	24.40	66.57	47.7

d_{xx} (μm) indicates the size of particle below which XX% of the sample lies (50% - median diameter)


Aldrich (310700):



d_{10} : 1.475 μm d_{50} : 4.302 μm d_{90} : 11.16 μm
<75% 7.022 μm <95% 14.53 μm
<10 μm 87.5%




Aldrich (243442):




d_{10} : 77.66 μm
<75% 156.8 μm
<10 μm 1.57%


d_{50} : 138.4 μm
<95% 183.2 μm





Aldrich (217646):



d_{10} : 0.577 μm d_{50} : 4.240 μm d_{90} : 16.04 μm
<75% 9.883 μm <95% 19.69 μm
<10 μm 75.3%




Acros (213490010):




10x



Wako (138-09675):




d_{10} : 2.743 μm d_{50} : 25.70 μm d_{90} : 63.10 μm
<75% 44.77 μm <95% 73.05 μm
<10 μm 21.6%



MnO₂ Old Sample:



10x



40x